

Evaluation of anti-CENP reactivity in samples exhibiting the centromere HEp-2 pattern, which is associated with a better prognosis within the limited cutaneous systemic sclerosis spectrum

Poster

Autoimmunity and Transplantation

Gerson D. Keppeke^{1, 2}, Diana Landoni^{2, 3}, Cristiane Kayser², Pedro Matos², Larissa Diogenes², Jessica Keppeke⁴, Silvia Helena Rodrigues², Luis Eduardo C. Andrade^{2, 4}.

(1) Departamento de Ciencias Biomédicas, Facultad de Medicina, Universidad Católica del Norte, Av Larrondo 1281, Coquimbo, Chile.

(2) Disciplina de Reumatologia, Departamento de Medicina, Universidade Federal de São Paulo, São Paulo, Brasil.

(3) Escuela de Graduados, Facultad de Medicina, Universidad de la República, Montevideo, Uruguay.

(4) Immunology Division, Fleury Medicine and Health Laboratories, Fleury group, São Paulo, Brasil.

gerson.keppeke@ucn.cl

Introduction/Objectives:

Anti-centromere antibodies (ACA) are among the most specific markers for systemic sclerosis (SSc), particularly the limited cutaneous form (lcSSc), where they are traditionally associated with a more favorable prognosis. The discrete speckled centromere pattern (AC-3) on HEp-2 cells strongly suggests the presence of antibodies directed against centromere proteins, especially CENP-B and CENP-A. However, the precise correlation of AC-3 reactivity with specific CENP antigens, as well as its clinical implications within a homogeneous lcSSc cohort, remains incompletely understood. This study aimed to investigate the frequency and distribution of anti-CENP-A and anti-CENP-B reactivities in sera from lcSSc patients displaying the AC-3 pattern, and to analyze their clinical and demographic associations, particularly in relation to organ involvement.

Methods:

A total of 76 patients fulfilling ACR/EULAR 2013 classification criteria for lcSSc were included. Among them, 38 exhibited the AC-3 pattern on HEp-2 IFA, while 48 were classified as Non-AC-3. Reactivity to CENP-B and CENP-A was evaluated by two solid-phase immunoassays: ELISA (for CENP-B) and line-blot (for both CENP-A and CENP-B). In parallel, clinical data from 68 lcSSc, including age, sex, disease duration, and major organ involvement (interstitial lung disease [ILD], pulmonary arterial hypertension, esophageal dysmotility, arthritis, and cardiac or renal manifestations) were reviewed. Statistical comparisons between groups and correlations between assays were performed.

Results:

Among the AC-3 samples, 92.1% were positive for CENP-B by ELISA and 84.2% by line-blot, while 81.6% reacted with CENP-A. Concordance for CENP-B reactivity between ELISA and line-blot was moderate (78.9%). Overall, 100% of AC-3-positive samples displayed reactivity to either CENP-A or CENP-B, confirming the high sensitivity of HEp-2 IFA in identifying ACA. Combined reactivity to both CENP-A and CENP-B in all platforms was observed in 68.4% of samples, whereas one serum exhibited exclusive anti-CENP-A reactivity, and seven were positive only for CENP-B. Clinically, AC-3-positive patients were older but showed significantly reduced frequency of ILD compared to Non-AC-3 lcSSc patients (10.5% vs. 54.2%; $p=0.001$). Frequencies of other organ manifestations, including esophageal involvement, arthritis, pulmonary hypertension, or renal and cardiac complications, were similar between groups.

Conclusions:

This study demonstrates that in lcSSc patients, the AC-3 centromere pattern on HEp-2 IFA is strongly predictive of anti-CENP reactivity, with CENP-B being the predominant target but anti-CENP-A also present, including rare cases of exclusive reactivity. Importantly, the presence of AC-3 in this lcSSc cohort was associated with a distinctly lower prevalence of ILD, reinforcing its role as a marker of milder disease course and more favorable prognosis. The incomplete concordance between ELISA and line-blot highlights the variability between solid-phase assays, underscoring the importance of HEp-2 IFA as a primary screening tool. Altogether, these findings support the clinical utility of combining immunofluorescence with antigen-specific assays for the precise characterization of anti-centromere reactivity in systemic sclerosis, and provide evidence that ANA patterns can further stratify lcSSc patients into a subgroup with reduced risk of severe pulmonary involvement.


key-points

- In a cohort composed exclusively of limited cutaneous Systemic Sclerosis patients, those with the AC-3 centromere HEp-2 pattern showed significantly lower prevalence of interstitial lung disease compared to those without AC-3.

- CENP-B is the dominant autoantigen in AC-3 samples, but anti-CENP-A reactivity was also prevalent and exclusive anti-CENP-A reactivity was observed in one sample.
- HEp-2 IFA was 100% sensitive for detecting anti-CENP antibodies, meaning all AC-3-positive samples reacted with CENP-B and/or CENP-A in at least one solid-phase immunoassay, reinforcing HEp-2 IFA as a highly reliable screening method.

Palabras Clave

(1) Systemic Scleroderma, (2) Autoantibodies, (3) Fluorescent Antibody Technique, (4) HEp-2 Cells, (5) Centromere.



Venn diagram for anti-CENP-B/A reactivity in ELISA and line-blot assays in samples from the AC-3 group.

Financiamiento:

FAPESP (Sao Paulo State Research Foundation) grant numbers #2017/20745-1,

#2021/04588-9 and #2023/17946-6, granted to GDK, LD and CK, respectively.

ANID #86220018, proyecto de Fortalecimiento del Doctorado en Ciencias Biomédicas - UCN.